

TismoAI Diagnostic Assistant for Instruments

Presented at Pittcon 2025

Tismo Technology Solutions (P) Ltd.

Bangalore, INDIA

Abstract

The rise of AI has revolutionized troubleshooting in industries handling complex instruments. Legacy systems evolved over decades, incorporating feature changes, bug fixes, and customer-driven updates, making efficient reference crucial. Modern instruments include detailed service manuals, firmware updates, and FAQs, yet troubleshooting remains highly complex.

This paper introduces an AI assistant that enhances accuracy and efficiency by diagnosing and resolving issues using diverse data sources and formats. With voice interaction capabilities, it empowers service engineers and customers, streamlining troubleshooting for analytical instruments. This AI-driven approach transforms maintenance workflows, reducing downtime and improving overall operational efficiency.

1. Introduction

Technical and analytical instruments play a critical role in medical, manufacturing, and communication technology. These sophisticated instruments require precise operation and maintenance, as failures can disrupt essential procedures. Troubleshooting them is often complex, time-consuming, and demands deep problem understanding. Service engineers must analyze vast knowledge from manuals, service reports, and past troubleshooting cases across various locations. They also face challenges in locating accurate information, deciphering random error codes, and piecing together system history to diagnose issues effectively.

The rapid advancement of AI technology presents an opportunity to transform troubleshooting. Generative AI enables applications to analyze vast data, perform complex analyses, and provide instant recommendations. This white paper outlines our experience in developing and deploying an AI-driven assistant for analytical instrument troubleshooting.

The AI assistant is designed to serve primarily two audiences—service engineers with technical knowledge and end-users who require step-by-step assistance. By

leveraging a comprehensive knowledge base, the assistant aims to bridge the gap between user queries and accurate solutions. Furthermore, the proposed AI assistant supports secure and flexible deployment options, ensuring usability across different operational environments.

This initiative not only addresses the inefficiencies of traditional troubleshooting methods but also sets a new standard for reliability, flexibility, and user experience in technical support systems.

2. Challenges of Manual

Troubleshooting

Manual troubleshooting of analytical instruments poses several challenges, including complexity of issues, human error, and dependence on expert knowledge. The process is often slow, error-prone, and costly, with limited access to critical data, further complicating resolution. This study analyzes these challenges to explore potential improvements.

2.1 Time-intensive Process

Technicians rely on systematic diagnosis and trial-and-error methods, making troubleshooting slow, especially for complex instruments. Identifying root causes, testing solutions, and addressing interdependencies take considerable time. This delays operations, reduces productivity, and impacts industrial throughput.

2.2 Dependency on Expert Knowledge

Sophisticated instruments require specialized knowledge, but experts are often unavailable, creating bottlenecks. Less experienced staff risk misdiagnosing issues, complicating repairs. The reliance on experts also makes scaling operations and managing instruments across locations difficult, highlighting the need for standardized troubleshooting approaches.

2.3 High Error Rates

Human errors are inevitable in manual troubleshooting, especially when interpreting complex data. Mistakes in fault identification, test analysis, or corrective actions prolong downtime and may damage instruments, increasing repair costs. Errors can also affect calibration, leading to inaccurate readings and unreliable analysis.

2.4 Limited Data Accessibility

updated repository is accessible through the assistant, allowing technicians to instantly retrieve accurate and comprehensive information. The assistant's robust search capabilities enable users to quickly find the required data, reducing inefficiencies from fragmented or outdated information and improving troubleshooting accuracy and speed.

3.2 24x7 Availability

Majority of the AI assistants are cloud-hosted and operate continuously. This feature can be leveraged to provide troubleshooting support anytime, anywhere. It




Figure 1: Key strategies for enhancing the troubleshooting experience

Effective troubleshooting requires access to historical and real-time instrument data. However, manual methods often involve fragmented systems and lack streamlined interfaces, making data retrieval difficult. Incomplete or delayed access hampers diagnostics and prevents root cause analysis. Without centralized data repositories, identifying recurring issues and implementing predictive maintenance becomes challenging.

3. Key Strategies

To overcome the challenges mentioned above, Tismo identified five key strategies to enhance the efficiency and effectiveness of troubleshooting analytical instruments by utilizing the capabilities of AI. These improvements can be achieved through the following implementations.

3.1 Centralized Data Source

AI assistants generally use a centralized repository to store essential instrument data, including manuals, error codes, and maintenance logs. This continuously

reduces the need for human intervention during off-hours, lowering operational costs while ensuring consistent service quality. By diagnosing issues instantly, it minimizes downtime and improves efficiency. Its constant availability delivers quick solutions, ensuring a reliable and stress-free troubleshooting experience, even in urgent or remote situations.

3.3 Voice Interaction

Modern browsers include built-in voice interaction support. When users grant microphone access, browsers utilize voice recognition to capture and process sound. With AI-powered assistance, hands-free troubleshooting becomes possible in challenging environments. This feature is especially useful for technicians in busy or confined spaces. The assistant interprets queries, provides step-by-step guidance, and confirms actions through spoken responses, reducing dependence on manual input. This voice interface enhances accessibility and efficiency, enabling seamless multitasking during troubleshooting.

3.4 Simplified Scalability

AI systems are designed to adapt effortlessly to growing data volumes and new products. Their architecture enables seamless integration of additional data sources, allowing them to evolve with organizational needs without added complexity or overhead. With access to extensive knowledge repositories, AI-powered assistants can operate at scale, supporting thousands of users simultaneously while maintaining optimal performance and response speed.

3.5 Consistency and Accuracy

By leveraging standardized processes and advanced AI algorithms, AI-powered systems eliminate human error, ensuring precise and consistent troubleshooting for complex instruments. They systematically analyze issues, process intricate queries, and retrieve the most relevant information, enabling faster and more accurate diagnostics. AI-driven approaches minimize variability in troubleshooting, reducing the risk of misinterpretation and ensuring reliable outcomes. As these systems continuously expand their knowledge base with newly added data, they improve over time, enhancing efficiency and effectiveness. This keeps the Assistant updated with latest data. This reliability fosters user confidence, as AI provides consistent, data-driven solutions, ultimately optimizing performance and minimizing downtime in critical applications.

4. Benefits

To achieve the key objectives, Tismo designed an AI assistant, the TismoAI Diagnostic Assistant or TADA. Tismo leveraged the key components and capabilities of AI, while also providing additional functionalities and tackling other potential shortcomings and limitations that are inherent in manual troubleshooting.

4.1 Unified Repository with Diverse data

TismoAI Diagnostic Assistant (TADA) brings together diverse data sources into a single framework, ensuring no critical information is overlooked. Its integration capabilities span across structured and unstructured formats, making it a powerful tool for end-to-end troubleshooting. Below mentioned are the different sources of information that are accepted by TismoAI Diagnostic Assistant, with some requiring additional preprocessing before uploading them to the repository to achieve the best query responses:

● Service Manuals and Technical Documentation:

Service manuals for analytical instruments are comprehensive documents provided by manufacturers,

outlining technical specifications, calibration procedures, maintenance protocols, and troubleshooting guidelines. They are essential for troubleshooting as they offer detailed, step-by-step diagnostics and instrument-specific solutions, ensuring accurate problem resolution and minimizing both downtime and operational errors.

● Issue Tracking Reports:

Issue tracking reports, such as those in Mantis, document software and hardware problems, including error logs, user-reported issues, and resolutions. They aid in troubleshooting analytical instruments by offering a historical record of recurring faults and successful fixes. This enables faster, data-driven diagnostics and reduces repetitive problem-solving efforts.

● Internal Knowledge Base Articles:

Internal knowledge base articles differ from service manuals as they are organization-specific and regularly updated. They provide real-world solutions, undocumented fixes, and insights from field engineers, making them more adaptable to issues not covered in official manuals.

● Email Threads:

Email threads aid troubleshooting by preserving real-time problem-solving discussions between service engineers and users. They offer direct insights into issues, field-tested solutions, and unique challenges, providing context-specific resolutions. Unlike static documentation, these exchanges capture evolving diagnostics, enabling faster, experience-driven responses. Relevant information from email threads can later be used to create internal knowledge base articles.

● Product Change Notes and Firmware Updates:

Product change notes (PCNs) outline modifications or improvements to a product, including hardware or software updates. Firmware updates involve changes to the firmware controlling the device. Both are essential for troubleshooting analytical instruments, as they offer insights into new features, bug fixes, and compatibility improvements, ensuring optimal performance.

● FAQs and Troubleshooting Guides:

FAQs and troubleshooting guides offer step-by-step solutions to common issues, helping users diagnose and resolve problems quickly. For analytical instruments, they provide predefined solutions, error explanations, and preventive measures, reducing downtime and minimizing expert intervention. This improves user efficiency and ensures accurate instrument operation.

FAQs help quickly find solutions to both recent and

relevance, ensuring users receive precise and actionable




Figure 2: Diverse sources of information supported by TismoAI Diagnostic Assistant

recurring issues. These documents are typically stored in CSV format.

● Solution Reports:

Solution reports document detailed troubleshooting steps, diagnoses, and resolutions for issues with an analytical instrument after production. They are ideal for troubleshooting, offering a systematic approach, insights into recurring issues, and helping technicians identify patterns. This reduces future troubleshooting time and ensures consistent instrument performance.

● Error Codes and Descriptions:

Error codes are essential for diagnosing and troubleshooting analytical instruments, providing precise information about malfunctions. They help users quickly identify issues, determine root causes, and apply corrective actions. By referencing error codes in manuals or support documentation, technicians can resolve problems efficiently, minimizing downtime and ensuring accurate operation.

4.2 Application of Generative AI

Generative AI in TismoAI Diagnostic Assistant creates tailored responses, guides, and solutions based on user queries' context. This personalized approach enhances

information. The assistant's responses have a conversational, human touch, improving the overall user experience.

4.3 LLM-Independent Architecture

TismoAI Diagnostic Assistant integrates seamlessly with multiple large language models (LLMs), offering organizations flexibility in selecting the model that best suits their needs. Whether prioritizing cost, security, performance, or compliance, businesses can configure TismoAI Diagnostic Assistant with the most appropriate LLM. For scalability and AI advancements, Google Gemini is ideal, while Amazon SageMaker is suitable for cloud-based, customizable solutions. Microsoft LUIS meets enterprise security and compliance needs, and OpenAI's GPT models excel in natural language fluency and innovation. This adaptability ensures cost efficiency, future-proofing, and alignment with organizational policies. TismoAI Diagnostic Assistant's architecture allows businesses to design analytical instruments to tailor AI support, optimizing troubleshooting efficiency and operational security.

4.4 Legacy System Support

TismoAI Diagnostic Assistant primarily supports modern instruments with internet connectivity, enabling real-time diagnostics and updates through IoT capabilities. These systems provide endpoints that allow the assistant to continuously monitor key parameters such as alarm status, temperature, and pressure, ensuring timely troubleshooting and proactive support.

For legacy instruments, the assistant integrates with adapters that connect these devices to online platforms. These adapters relay data, including instrument status, alarms, temperature, and pressure readings, to the AI assistant via various communication protocols. This ensures both modern and legacy instruments benefit from AI-driven diagnostics and real-time monitoring, regardless of their age or design.

4.5 Frequently Asked Questions (FAQs)

– One Click Away

TismoAI Diagnostic Assistant features a robust FAQ system that delivers quick answers to common problems, allowing users to resolve issues independently. This reduces the burden on support resources, freeing them up for more complex cases. With a single click, the assistant provides a detailed solution, significantly reducing troubleshooting time compared to traditional methods.

4.6 Intuitive User Experience

TismoAI Diagnostic Assistant offers an intuitive interface and contextual understanding, simplifying troubleshooting for users of all expertise levels. Its conversational and interactive features create a seamless, engaging experience, enhancing user satisfaction. A chatbot popup will be integrated on the analytical instrument website, allowing users to quickly query a frequently asked question or troubleshoot a new issue. The assistant's responsive UI ensures smooth functionality on both desktops and smartphones.

4.7 Access to Previous Sessions

Each user session with the assistant generates a unique public session ID, enabling users to easily return to previous conversations. A concise session summary is automatically created and used as context during revisit, allowing the assistant to quickly understand past interactions. This eliminates the need for users to repeat details, streamlining the support process. Especially when troubleshooting with a service engineer, TismoAI Diagnostic Assistant saves valuable time and improves

issue resolution, making the experience more efficient and user-friendly.

4.8 Reduced Downtime

TismoAI Diagnostic Assistant's scalability ensures immediate, effective solutions, greatly reducing the time needed to diagnose and resolve issues. Unlike traditional support systems that require waiting in queues or contacting customer support within specific time zones, TismoAI Diagnostic Assistant offers real-time, 24/7 assistance. This eliminates delays, enabling users to troubleshoot problems instantly, whether in a local lab or globally. By removing these bottlenecks, businesses experience minimal downtime, ensuring seamless operations. This enhances productivity and maximizes profitability by keeping analytical instruments running efficiently without prolonged disruptions.

5. Implementation

The implementation of TismoAI Diagnostic Assistant is centered on a secure yet decoupled architecture, ensuring seamless integration with a large language model (LLM) while maintaining system modularity. This architectural design enables the deployed LLM to retrieve relevant insights from an internal repository, which houses extensive domain-specific data related to the analytical instrument. The repository accommodates both structured and unstructured data, necessitating rigorous preprocessing techniques to refine and filter the dataset. This preprocessing ensures that only the most relevant information is retained, optimizing the quality and accuracy of query responses. Additionally, a dedicated alarm server continuously monitors the real-time status of the analytical instrument, providing live updates on operational conditions. All processed data and real-time alerts are seamlessly presented through an AI Assistant User Interface—a smart chat assistant deployed on the website—designed to facilitate intuitive and efficient interaction with users. The following sections provide a detailed analysis of each architectural component and its role in enhancing system functionality:

5.1 Data Processing

Data preprocessing is a crucial step in refining raw textual information into a structured and meaningful format, enabling efficient retrieval and analysis. When dealing with large volumes of technical documentation—ranging from service manuals to issue tracking reports and email threads—raw data often contains inconsistencies, redundancies, and irrelevant information that can hinder accurate AI-based

responses. Preprocessing involves techniques such as text extraction, segmentation, standardization, entity recognition, and summarization to ensure that the data is clean, structured, and contextually relevant. Without preprocessing, an AI system may struggle with retrieving precise answers, leading to ambiguous or irrelevant responses. By implementing a well-defined preprocessing pipeline, the data fed into the internal repository becomes highly optimized, significantly improving the performance of the deployed large language model (LLM).

To structure the data effectively, different

signatures and disclaimers, and topic modeling to categorize conversations based on subject matter.

Documents related to product updates and troubleshooting require additional refinement. Product change notes and firmware updates undergo version control tracking, ensuring that modifications are linked to specific firmware versions, while dependency mapping highlights affected components. FAQs and troubleshooting guides benefit from question-answer pair extraction, structuring content into easily retrievable sections for AI-driven responses. Solution reports, which document diagnostic findings and




Figure 3: TismoAI Diagnostic Assistant Architecture Diagram

preprocessing techniques are applied based on document type. Service manuals, for example, undergo OCR processing if they are in scanned PDF format, followed by text segmentation to extract key sections such as installation guidelines, maintenance procedures, and troubleshooting steps. Issue tracking reports, which document past failures and solutions, require deduplication to remove redundant records, standardization for consistent terminology, and sentiment analysis to categorize issues by severity. Meanwhile, internal knowledge base articles benefit from content normalization and entity linking, connecting relevant information across different sources. Email threads, which often contain fragmented discussions, undergo thread collapsing to group messages by topic, noise removal to eliminate

resolutions, require pattern recognition to identify recurring issues, root cause analysis extraction to highlight key findings, and cross-referencing to link reports with related technical documentation. These preprocessing steps collectively enhance the accessibility and usability of data, ensuring that AI-generated responses are both accurate and contextually relevant.

The impact of preprocessing can be illustrated with a simple example. Suppose an engineer queries the AI assistant about a recurring calibration issue. Raw data, such as an unstructured issue tracking report, might contain multiple instances of the same problem, different terminologies used by various engineers, and lengthy email threads with scattered information. This could lead the LLM to generate an unclear or

incomplete response. In contrast, preprocessed data would have a standardized issue description, structured failure reports, and consolidated solution records, enabling the LLM to quickly retrieve the most relevant and precise response. As a result, preprocessing ensures that AI-driven troubleshooting is efficient, reducing ambiguity and improving problem resolution times for users of the analytical instrument.

5.2 Internal Repository

TismoAI Diagnostic Assistant evaluated two different approaches for managing its internal repository: a data store managed by Google and a vector database entirely controlled by Tismo. The first approach leverages Google’s internal data storage system within Vertex AI,

they efficiently store and retrieve knowledge, allowing AI models to generate highly relevant responses based on user queries.

Pinecone is a leading vector database that enables real-time similarity searches by indexing vector embeddings. These embeddings represent the underlying meaning of data rather than exact keyword matches, allowing for more accurate and contextually relevant responses. When a user submits a query, Pinecone searches for the closest vector representations, ensuring precise retrieval of information. This capability is crucial in AI-powered troubleshooting systems like TismoAI Diagnostic Assistant, where accurate knowledge retrieval significantly enhances response quality.




Figure 4: Indexing in Vector Database

while the second employs Pinecone, a vector database optimized for AI-driven applications.

A vector database (Vector DB) is a specialized storage system designed for managing high-dimensional vector embeddings rather than traditional tabular data. Unlike conventional relational databases that rely on structured query language (SQL) for retrieving information, vector databases enable semantic search by storing and querying data in vectorized form. This makes them ideal for AI applications, especially when handling unstructured data such as text, images, and audio. Vector DBs are preferred for LLM training because

In contrast, Google’s Vertex AI leverages its own internal data store, incorporating advanced AI-driven indexing and processing techniques. Google performs extensive pre-processing, improving the efficiency of knowledge retrieval by integrating structured and unstructured data seamlessly. One of Google’s key advantages over Pinecone is its hybrid search mechanism, which combines keyword-based search with vector-based semantic search. This hybrid approach improves retrieval accuracy by balancing structured query processing with deep semantic understanding. However, Google’s approach is less configurable compared to Pinecone, as it abstracts



Figure 5: LLM Agnostic Design of TismoAI Diagnostic Assistant

much of the indexing and optimization behind proprietary algorithms.

While Pinecone’s reliance on pure semantic search allows for a more tailored and fine-tuned retrieval mechanism, Google’s hybrid search method enhances response accuracy for a broader range of queries. Organizations leveraging TismoAI Diagnostic Assistant can choose between these approaches based on their need for configurability, response precision, and AI workflow integration.

5.3 LLM Environment

To demonstrate the LLM-agnostic nature of our AI assistant, two distinct approaches were explored by Tismo: integration with OpenAI’s GPT 4o-mini alongside a vector database, and Google Vertex AI combined with Dialogflow. These implementations highlight the adaptability of TismoAI Diagnostic Assistant across different ecosystems. This also ensures flexibility, cost optimization and regulatory compliance.

- **Integration with OpenAI GPT 4o-mini and Pinecone Vector Database:**

OpenAI’s GPT-4o mini, when coupled with a vector database such as Pinecone, enhances retrieval-augmented generation (RAG) capabilities. This integration enables efficient query processing by leveraging Pinecone’s semantic search over extensive document repositories, ensuring highly relevant responses. GPT-4o Mini further improves real-time AI assistance with low-latency query processing, allowing for instant and accurate information retrieval. Additionally, the expandable knowledge base powered by vector embeddings allows continuous updates to the AI’s knowledge repository without requiring full retraining.

- **Integration with Google Vertex AI and Dialogflow**

TismoAI Diagnostic Assistant was also integrated with Google’s Vertex AI and Dialogflow to leverage their robust cloud-based AI capabilities. This integration provides scalability and model customization by enabling organizations to fine-tune and deploy custom LLMs tailored to their specific use cases. Seamless integration with enterprise cloud infrastructure ensures high availability and data security. Additionally,

Dialogflow enhances AI-driven conversational experiences, offering advanced assistant functionalities.

5.4 Alarm Server

TismoAI Diagnostic Assistant is designed to seamlessly integrate with both modern and legacy analytical instruments, enabling real-time monitoring and alarm status retrieval without necessitating modifications to the instrument's configuration. For modern instruments equipped with internet connectivity, TismoAI Diagnostic Assistant leverages pre-built API endpoints to query the alarm status instantaneously. This capability allows users to monitor the operational health of their instruments within a laboratory setting or across global facilities. Users can directly query the instrument's alarm status through the TismoAI Diagnostic Assistant, either verbally or via text input, enabling quick access to critical information without navigating complex interfaces. To ensure data security and compliance, TismoAI Diagnostic Assistant strictly adheres to established security protocols by querying only the necessary data while discarding any extraneous information. Additionally, the system is designed to comply with industry regulations such as FDA 21 CFR Part 11, and Good Manufacturing Practice (GMP) guidelines, ensuring data integrity, traceability, and audit readiness. This approach not only facilitates efficient troubleshooting but also enhances predictive maintenance by providing real-time insights into instrument status.

For legacy instruments that lack direct internet connectivity, TismoAI Diagnostic Assistant employs a specialized adapter tailored to the instrument's configuration and accessibility constraints. This adapter selects an appropriate communication protocol that aligns with the manufacturer's specifications and the organization's security policies, ensuring seamless integration. The adapter can be designed to interface exclusively with a specific instrument model or to support a broader range of instruments, providing flexibility for organizations in deploying their monitoring solutions. Additionally, access control mechanisms enable organizations to restrict alarm status queries to authorized individuals or laboratories, thereby maintaining data integrity and security. As with modern instruments, real-time alarm status retrieval from legacy systems enhances troubleshooting efficiency and supports data-driven decision-making in laboratory management.

5.5 AI Assistant User Interface

TismoAI Diagnostic Assistant is designed to function as a chatbot popup, enabling effortless integration into existing websites without necessitating extensive code

modifications. The user interface (UI) is highly responsive and intuitive, ensuring a seamless experience across various devices, including smartphones, tablets, and desktops. A proprietary TismoAI Diagnostic Assistant icon is positioned at the bottom-right corner of the screen, providing quick access to the assistant. The UI is both expandable and collapsible, allowing users to minimize or close the assistant as needed.

To enhance usability, the assistant initially displays four frequently asked questions (FAQs) in button format at the top of the popup. Users can select any of these options to receive instant, curated responses. Additionally, upon activation, the assistant inquires whether the user would like to resume a previous conversation. If the user opts to continue an earlier session, they are prompted to enter a session ID corresponding to the desired conversation. Upon verification, the AI retrieves the internally stored summary of the previous session, establishing contextual continuity and enabling a seamless transition.

If the user prefers to initiate a new conversation, they may indicate so by selecting the appropriate option. The assistant supports both textual and voice-based interactions, allowing users to pose queries, ask follow-up questions, or introduce entirely new topics. The AI is capable of discerning contextual shifts, ensuring relevant and coherent responses to diverse inquiries.

Upon concluding their interaction, users have the option to either close the popup, terminating the session, or issue a "summarize" command. This action generates a summary of the current conversation and assigns a unique session ID, which can be utilized for future reference. Notably, minimizing the popup retains the conversation in an active state, whereas closing it completely ends the session. This approach ensures that users can flexibly engage with the assistant while maintaining control over their interactions.

6. Upcoming Capabilities

TismoAI Diagnostic Assistant is set for an exciting future, continuously evolving to incorporate more advanced features. As technology progresses, it will adapt to meet growing user demands, improving efficiency, accuracy, and performance. With ongoing enhancements to LLM models, the assistant will continue to improve over time. Its flexible architecture ensures seamless integration, making it a dynamic, future-proof system. Beyond the features discussed, there is potential to develop additional capabilities tailored to an organization's needs, helping the assistant integrate into their existing practices. A few of these are outlined below:

6.1 On-Premise Deployment

TismoAI Diagnostic Assistant will enable local deployment on analytical instruments, providing real-time AI assistance without the need for constant internet access. By running an LLM on organization servers, users will be able to access troubleshooting support even in offline or restricted environments. This solution will ensure faster response times, improved data security, and uninterrupted assistance. Local processing will reduce reliance on cloud services, making it an ideal choice for industries focused on data privacy and reliability, especially in remote or high-security settings.

6.2 Secure Email-Based Session Management

TismoAI Diagnostic Assistant will introduce email-based session management, adding an extra layer of security and accessibility to past interactions. Users will receive a session code via email, allowing them to securely revisit previous troubleshooting conversations without losing context. This enhancement will protect user data, prevent unauthorized access, and streamline issue resolution by maintaining a record of past discussions. By implementing this feature, TismoAI

6.3 Smart Feedback Loop

A comprehensive feedback mechanism will be introduced to refine TismoAI Diagnostic Assistant's responses based on user interactions. Users will be able to rate responses, flag inaccuracies, and suggest improvements, fostering a dynamic learning loop. This feedback will be analyzed to enhance response accuracy, ensuring the assistant continuously evolves to meet user needs. By incorporating machine learning-driven insights, TismoAI Diagnostic Assistant will become more precise over time, reducing errors and providing increasingly relevant solutions. This will ultimately boost user confidence and satisfaction, ensuring a more effective troubleshooting experience.

6.4 Multilingual Text and Audio Support

To serve a global audience, TismoAI Diagnostic Assistant will feature expanded multi-language support, offering both textual and audio-based responses. Users will receive real-time assistance in their preferred language, eliminating language barriers and improving accessibility. Advanced speech synthesis and recognition will enable seamless voice interactions,




Figure 6: TismoAI Diagnostic Assistant in action

Diagnostic Assistant will ensure that users can securely track and continue their troubleshooting journey without starting from scratch, enhancing both security and convenience.

making troubleshooting more intuitive. This enhancement will benefit multilingual teams, reducing miscommunication and ensuring a smoother troubleshooting experience across diverse regions.

Ultimately, it will boost adoption and efficiency in international markets, broadening TismoAI Diagnostic Assistant's reach and impact.

6.5 CDN-Powered Assistant for Seamless Web Deployment

TismoAI Diagnostic Assistant will integrate with a Content Delivery Network (CDN), enabling the assistant to run efficiently on any website worldwide. This integration will ensure low-latency performance, faster response times, and seamless accessibility across multiple geographies. By leveraging CDN technology, businesses will be able to embed TismoAI Diagnostic Assistant into their websites, providing instant AI-driven assistance to users without complex installations. This scalability will make it an ideal solution for customer support, product troubleshooting, and knowledge-sharing platforms, ensuring a consistent and high-performance experience for users across the globe.

7. Conclusion

The emergence of AI-based troubleshooting software, as represented by TismoAI Diagnostic Assistant, represents a turning point in analytical instrument maintenance and support. By resolving key issues such as the spread of knowledge, human error, and reliance on expert technicians, TismoAI Diagnostic Assistant presents a new model of enhanced efficiency, precision, and accessibility. With its centralized knowledge base, real-time diagnosis, voice-based interface, and architecture independent of large language models, the AI assistant simplifies the troubleshooting process, minimizing downtime and enhancing continuity of operations.

Moreover, the system's flexibility ensures that new and existing instruments alike can leverage AI-based diagnostics, thereby closing the gap between conventional and contemporary maintenance procedures. With upcoming features such as on-premises deployment, multilingual support, and intelligent feedback mechanisms, TismoAI Diagnostic Assistant continues to innovate, constantly enhancing its functionalities to align with industry needs.

As the evolution of AI technology continues, its application in technical support and problem-solving is bound to become increasingly significant. By leveraging the potential of AI, organizations become well positioned not only to transform their troubleshooting processes but also to establish new standards for reliability, scalability, and user experience. The future of instrument maintenance is moving towards intelligent, data-driven solutions, with

TismoAI Diagnostic Assistant spearheading the revolution in this paradigm shift.

8. References

- [1] Kou Murayama, Adam Blake, Tyson Kerr and Alan Castel, "When Enough Is Not Enough", AMA. 2015. [Online]. Available: <https://psycnet.apa.org/record/2015-52761-001>
- [2] Yang Zhang, Yufei Wang and Lina Yao, "When Large Language Models Meet Vector Databases: A Survey", AMA. 2023. [Online]. Available: <https://arxiv.org/html/2402.01763v1>
- [3] Amir Elmishali, Roni Stern and Meir Kalech, "An Artificial Intelligence paradigm for troubleshooting software bugs", AMA. 2018. [Online]. Available: <https://doi.org/10.1016/j.engappai.2017.12.011>
- [4] Krishnaram Kenthapadi and Ankur Taly, "Grounding and Evaluation for Large Language Models: Practical Challenges and Lessons Learned (Survey)", ACM. 2024. [Online]. Available: <https://dl.acm.org/doi/10.1145/3637528.3671467>